

Delivering nasal high flow with Optiflow™ at home  
and in long-term care




# Introducing the F&P myAIRVO™ 2, delivering respiratory support and airway hydration with a focus on comfort and ease of use.

The F&P myAIRVO 2 has been designed specifically to be an integrated solution for delivering nasal high flow with Optiflow™. Fisher & Paykel Healthcare's leading humidification technology allows the AIRVO 2 to comfortably provide high flows of air/oxygen mixtures to spontaneously breathing patients, through the unique Optiflow nasal interface.


## Optiflow™+ Nasal Interface

- Soft, flexible prongs
- Wide-bore design reduces gas jetting
- Prongs contoured to the patient's nose



## Adjustable temperature and flow settings

- Three temperature settings (37, 34, 31 °C) help achieve comfort and compliance
- Integrated flow generator delivers a wide flow range (2 - 60 L/min) - no wall air supply required



## Heated breathing tube

- Dual spiral heater wires and unique integrated temperature sensor
- No separate temperature probes or heater-wire adapters required

## Supplementary oxygen (when required)

- Oxygen can be added from a wall supply or cylinder
- Inbuilt ultrasonic oxygen analyzer requires no calibration, service or replacement

Product codes	
PT100US	myAIRVO 2 humidifier
OPT942E	Optiflow+ small interface (1-pack)
OPT944E	Optiflow+ medium interface (1-pack)
OPT946E	Optiflow+ large interface (1-pack)
OPT970E	Optiflow+ Tracheostomy interface (1-pack)
OPT980E	Optiflow+ mask interface adapter (1-pack)
900PT500E	Heated breathing tube (1-pack)
900PT500	Heated breathing tube (10-pack)
900PT290E	MR290 chamber (1-pack)
HC360	Reusable water chamber
900PT400	Compact stand
900PT401	Water bag (2-pack)

Find out more at [www.fphcare.com/myairvo](http://www.fphcare.com/myairvo)

**Fisher & Paykel**  
HEALTHCARE